

Instructions for
Colloidal Staining with AzurGel-K

for Gel format 10 x 10 x 0,1 cm

Cat No.: GF 10002

- **Basics**

This staining procedure uses the colloidal properties of Coomassie™ Brilliant Blue G 250 dyes in the presence of inorganic acid and a high concentration of salt. Optimal staining of proteins on a clear background can be achieved with this kit.

- **References**

Neuhoff V, Arold N, Taube D, Ehrhardt W (1988) *Electrophoresis* 9, 255 - 262

Neuhoff V, Stamm R, Eibl, H (1985) *Electrophoresis* 6, 427 - 448

- **Advantages**

Higher sensitivity as customary Coomassie™ Brilliant Blue dyes

Protein staining through the whole cross-cut. (< 1 mm)

Waterclear background after destaining

Stable fixing after staining

Suitable for Tris/Glycine/SDS- and Tris/Tricine-gels

Fewer changes of solution

Solvent free destaining using D.I.-water

- **Specification**

Stainer A: 500 ml, contains (Ammonium sulfate and Phosphoric acid)

Stainer B : 125 ml, contains Coomassie™ Brilliant Blue G 250

which is enough for 25 gels, 10 x 10 cm

Storage: room temperature

Shelf Life: 6 month, see package

Detection limit: 6 ng BSA (with 5 mm well widths and 1 mm gel diameter)

Duaration: ~ 12 h, 3 – 4 steps

Handling: Wear eye protector and gloves while handling the gels and solutions.

- **General Guidelines**

- **Shaker**

Accomplish fixing and staining with a shaker by 50- 70 rounds per minute.

- **Staining Containers**

Only use clean, dense staining containers to avoid evaporation of methanol. Ensure that the container diameter and depth is sufficient to permit gel coverage with 100 ml of solution per gel. If possible, stain only one gel per staining dish, but not more than two gels per staining dish.

- **Tris/Glycine/SDS gels**

1. Staining

Before staining, prepare the following solution:

| | |
|-------|--|
| 55 ml | D.I. water |
| 20 ml | methanol |
| 20 ml | Stainer A |
| 5 ml | Stainer B (shake well before use) |

Incubate the gel on a shaker in the prepared solution for several hours or overnight. After 3 hours, the bands are visible, but the complete staining of the gel cross-cut is finished after 8.5 hours. The background of a gel with a low concentration of acrylamide (< 10 % T) is often more stained, because the colloidal particles enter into the large polyacrylamide pores of these gels.

2. Washing/Destaining

Put the gel in clean staining bowl and fill up to 200 ml with D.I. water. Shake the gel until the background of the gel becomes waterclear. Optimal destaining is reached after 7 hours, nevertheless it is possible to leave the gel in the water for up to three days (ie. over the weekend) without the loss of protein band intensity. Obstinate background staining can be cleaned with 250 ml/l methanol. Long term storage (more than 3 days) can cause a decrease of the protein bands intensity.

3. Storage

If you want to store the gel for a long time in a solution, add 200g/l Ammonium sulfate to the solution. This process eliminates the colour compounds from passing over to the solution and the bands remain stained.

4. Drying

If the gel is to be kept for documentation or other reasons, usage of AnaFix gel drying kit is recommended. The gels are dried by evaporation on the lab bench without shrinking. When staining with AzurGel-K time of incubation with AnaRapid gel drying solution should not exceed 5 –15 minutes or the destaining of the protein bands may occur.

- **Tris/Tricine-Gels**

- 1. Fixing**

Incubate the gel for 15 minutes in the following solution:

| | |
|---------|----------------------------|
| 40 ml/l | D.I. water |
| 50 ml/l | methanol |
| 10 ml/l | acetic acid (concentrated) |

- 2. Staining**

Incubate the gel after fixing for 15 minutes in the following solution:

| | |
|-------|------------|
| 55 ml | D.I. water |
| 20 ml | methanol |
| 20 ml | Stainer A |

Afterwards, add 5 ml Stainer B (**shake well before use**) and incubate the gel on a shaker for several hours or overnight in this solution. After 3 hours bands are visible, but the complete staining of the gel cross- cut is accomplished after 8.5 hours.

- 3. Washing/Destaining**

Put the gel in a second, clean staining bowl and wash it with 200 ml D.I. water. Then apply the same procedures as described in the Tris/Glycine instruction.

- **Cleaning**

10 ml HCl 32 %
90 ml methanol

This solution can be stored and may be used several times.

- **Trouble shooting**

| Problem | reason/cause | action taken |
|---|--|---|
| Staining solutions contains solid particles of stainer compound | Aggregation of colloidal particles | This is normal – no action required |
| Gel matrix completely stained blue | Not enough Stainer A used, staining compound therefore not in colloidal format | Destain gel in 25 % methanol, stain again if necessary |
| Proteins after 12 hours of staining show weak staining effect | Not enough Methanol or evaporated | Most of staining compound too strong in colloidal solution; staining time is extremely extended |
| Bands from different proteins with similar amounts show variations in staining effectiveness. | Different strongness in affinity of the stainer is an obviously known effect | For quantifications, always use calibrated amounts of proteins of interest |