

Vertikal Isoelectric Focussing

<i>Sample Buffer</i>	2-fold concentrate, for preparing 10 ml: 1,0 ml Ampholyte solution 3,0 ml Glycerol 5 mg Phenol Red 6,0 ml Water, ultra pure	
<i>Anoden Buffer I</i>	ready-to-use solution, also suitable for reverse IEF 0,34 % L-Asparatic acid 0,36 % L-Glutamic acid	
<i>Anode Buffer II</i>	50-fold concentrate 2,4 % Phosphoric acid, 85 %	
<i>Cathode Buffer I</i>	for IEF pH 3 – 10, ready-to-use solution 0,35 % Arginine base 0,29 % Lysine base	
<i>Cathode Buffer II</i>	for IEF 3 – 7 und other narrow ranges, ready-to-use solution 0,58 % Lysine base	
<i>Electrophoresis</i>	pH 3 – 10 and pH 3 – 7	pH 5 – 7 and 4 – 6 etc. 100 V, 60 min 100 V, 30 min 200 V, 60 min 200 V, 30 min 500 V, 30 min 500 V, 120 min
<i>Detection</i>	There are several methods for fixing and detection of protein bands, please find one below. Fixing: 20 % Trichloroacetic acid, 15 min, RT Rinsing: Water, ultra pure, 2 x 5 min Staining: 0,1 % Serva Violet 17 in 10 % Phosphoric acid, 15 min, RT Please prepare solution freshly. Destaining: Water, ultra pure, 3 x 10 min	