

## *Electrophoresis Protocol for anamed Tris Glycine Gels*

### *Native Electrophoresis*

#### Sample Buffer, 2-fold

Tris/Cl <sup>-</sup>	200	mM
Glycerol	20	%
Bromophenol Blue	0,005	%
pH	8,6	

#### Running Buffer, 10-fold

Tris	250	mM
Glycine	1,92	M
pH	8,5	

Do not use acid or base to adjust pH.

#### Running Conditions

Voltage	125	V	constant
Expected Current	~ 12	mA	per gel at the beginning of the run
	~ 6	mA	per gel at the end
Run Time	1 – 12	h	

### *Denaturing (SDS) Electrophoresis*

#### Sample Buffer, 2-fold (for reducing conditions reducing agent can be added)

Tris/Cl <sup>-</sup>	125	mM
Glycerol	20	%
SDS	4	%
Bromophenol Blue	0,005	%
pH	6,8	

#### Running Buffer, 10-fold

Tris	250	mM
Glycine	1,92	M
SDS	1	%
pH	8,3	

Do not use acid or base to adjust pH.

#### Running Conditions

Voltage	125	V	constant
Expected Current	~ 40	mA	per gel at the beginning of the run
	~ 20	mA	per gel at the end
Run Time	~ 90	min	